

ACCUPET



Revolutionary results come from evolutionary products



EVOLUTION



Ultra lightweight durable construction



Color coding for easy identification



Fully autoclavable

121° C



UV Resistant



Include ISO 17025-:2005 accredited calibration certificate



3 Years Warranty



Manifold rotates 360° to allow any convenient pipetting angle

Specially designed universal tip fitting nose cones for easy tip loading and ejection

The thermoplastic elastomer handle prevents transfer of body heat to the pipette to insure high accuracy and adds comfort to your grip even during continuous use.



- Low force pipetting
- Ergonomic design
- Light & smooth punger action
- Soft feel handle grip
- Suitable for both - left and right hand users

Soft touch tip ejection

Locking mechanism for preventing accidental volume changes

- Large and clear 4 digit display with smaller increments for wider volume selection options
- Continuously adjustable indicator
- Volume display is always visible while pipetting

Universal tip shaft

Revolutionary results come from evolutionary products



EVOLUTION

Variable Volume Micropipette

Cat. No.	Volume Range [μ l]	Increments [μ l]
EV-2.5	0.1-2.5	0.002
EV-10	0.5-10	0.02
EV-20	2-20	0.02
EV-50	5-50	0.1
EV-100	10-100	0.1
EV-200	20-200	0.2
EV-1000	100-1000	1
EV-5000	500-5000	10
EV-10ML	1000-10000	20

Fixed Volume Micropipettes

Cat. No.	Volume [μ l]	Cat. No.	Volume [μ l]
EF-2.5	2.5	EF-200	200
EF-5	5	EF-250	250
EF-10	10	EF-500	500
EF-20	20	EF-1000	1000
EF-25	25	EF-2000	2000
EF-50	50	EF-5000	5000
EF-100	100	EF-10ML	10000

8 Channel Micropipettes

Cat. No.	Volume Range [μ l]	Increments [μ l]
EV8-10	0.5 - 10 μ l	0.02
EV8-20	2 - 20 μ l	0.02
EV8-50	5 - 50 μ l	0.1
EV8-100	10 - 100 μ l	0.1
EV8-200	20 - 200 μ l	0.2
EV8-300	30 - 300 μ l	0.2

12 Channel Micropipettes

Cat. No.	Volume Range [μ l]	Increments [μ l]
EV12-10	0.5 - 10 μ l	0.02
EV12-20	2 - 20 μ l	0.02
EV12-50	5 - 50 μ l	0.1
EV12-100	10 - 100 μ l	0.1
EV12-200	20 - 200 μ l	0.2
EV12-300	30 - 300 μ l	0.2

Pipette Stand

Cat. No.	Description
SAS-01	Carousel stand for 6 pipettes
SAS-02	Shelf mount pipette holder

- Powerful yet quiet pump can fill 25 ml pipette under 4 seconds
- Soft grip push buttons
- Autoclavable nose cone
- Double Safety valve and hydrophobic filters provide double protection against fluid penetration
- Two different speed modes – High, Low along with Gravity drain.
- The pump speed can also be fine-tuned by varying finger pressure on operating knobs for better control of speed
- Cadmium Free environment friendly NiMH batteries
- Batteries can be changed very easily by the end user
- The intelligent charger prevents over charging / heating of batteries
- Low battery indicator
- The specially designed desktop stand enables the PipetHelp to be charged while resting on it
- The PipetHelp can hold the serological pipette while resting on the desktop stand



Pipette Controller

Cat. No.	Description
PH01-R	AccuHELP - Red
PH01-B	AccuHELP - Blue
PH01-O	AccuHELP - Orange
PH01-G	AccuHELP - Green

CAESA-LAB INC.

Tel: (450) 441-5408

www.caesalabinc.com

info@caesalabinc.com