

Collection & Containment Solutions



Autoclavable Biohazardous Waste Bags

Autoclavable Biohazardous Waste Bags are built to withstand the rigors of Steam Sterilization while meeting the requirements for the collection, storage, transport, and processing of biohazardous waste within the acute care facility.

Autoclavable bags made from high strength Polypropylene (PP) Film have an operating temperature of 274°F (134°C) and can tolerate temperatures of up to 285°F (140°C). High Density Polyethylene (HDPE) film autoclavable bags have an operating temperature of 250°F (120°C) and can tolerate temperatures of up to 260°F (127°C).

Autoclavable Biohazardous Waste Bags - HDPE Film w/ Indicator, Flat Pack, Twist Ties

Part #	Size/Width	Gauge	Color	Capacity	Qty/case
Print Label: Max Temp 260° F Biohazard Labeling; Biohazard Symbol					
<u>8-900</u>	8.5x11	25.4 mic	Red/Black	1-2 gal	1000
<u>8-904</u>	14x19	25.4 mic	Red/Black	3-4 gal	200
<u>8-908</u>	19x24	25.4 mic	Red/Black	7-10 gal	200
<u>8-912</u>	25x35	25.4 mic	Red/Black	13-16 gal	200

Autoclavable Biohazardous Waste Bags - Polypropylene Film w/ Indicator, Coreless Roll, Twist Ties

Part #	Size/Width	Gauge	Color	Capacity	Qty/case	Qty/pkg
Print Label: Biohazard Labeling; Biohazard Symbol						
<u>D5124</u>	25x35	1.75 mil	Buff/Black	13-16 gal	200	20

Autoclavable Biohazardous Waste Bags - Polypropylene Film w/ Indicator, Flat Pack, Twist Ties

Part #	Size/Width	Gauge	Color	Capacity	Qty/case
Print Label: Biohazard Labeling; Biohazard Symbol; English/French					
<u>5064.F</u>	30x38	2 mil	Yellow/Black	20-30 gal	125
<u>5142.F</u>	24x36	1.75 mil	Clear/Black	12-16 gal	200
<u>842Y</u>	19x24	2.2 mil	Yellow/Black	7-10 gal	200
<u>844Y</u>	25x30	2.2 mil	Yellow/Black	12-14 gal	200
<u>846Y</u>	25x35	2.2 mil	Yellow/Black	13-16 gal	200
<u>848Y</u>	38x48	2.2 mil	Yellow/Black	44 gal	100
<u>856Y</u>	31x38	2.2 mil	Yellow/Black	20-30 gal	100

Print Label: Max. Temp 285°F. Biohazard Labeling, Biohazard Symbol, English/Spanish



<u>5001.4</u>	8x12	2 mil	Red/Black	1-2 gal	400
<u>5010.1</u>	12x24	2 mil	Red/Black	2-3 gal	100
<u>5012</u>	14x19	2 mil	Red/Black	3-4 gal	500
<u>5013</u>	19x24	2 mil	Red/Black	7-10 gal	200
<u>5028</u>	40x70	2 mil	Red/Black	55-60 gal	100
<u>5119</u>	19x23	1.75 mil	Buff/Black	7 gal	200
<u>5237</u>	25x35	2 mil	Red/Black	13-16 gal	200
<u>8-300</u>	36x35	1.8 mil	Clear/Black	33 gal	100
<u>842</u>	19x24	1.8 mil	Red/Black	7-10 gal	200
<u>844</u>	25x30	1.8 mil	Red/Black	12-14 gal	200
<u>846</u>	25x35	1.8 mil	Red/Black	13-16 gal	200
<u>848</u>	38x47	1.8 mil	Red/Black	44 gal	100
<u>852</u>	12x24	1.8 mil	Clear/Black	2-3 gal	200
<u>854</u>	25x30	1.8 mil	Clear/Black	12-14 gal	200
<u>855</u>	31x38	1.8 mil	Clear/Black	20-30 gal	100
<u>856</u>	31x38	1.8 mil	REd/Black	20-30 gal	100
<u>858</u>	31x38	1.8 mil	Buff/Black	20-30 gal	100
<u>864</u>	38x47	1.8 mil	Buff/Black	44 gal	100
<u>872</u>	38x47	1.8 mil	Clear/Black	44 gal	100
<u>880</u>	8x12	1.8 mil	Orange/Black	1-2 gal	400
<u>881</u>	14x18	1.8 mil	Orange/Black	3-4 gal	400
<u>882</u>	19x24	1.8 mil	Orange/Black	7-10 gal	200
<u>884</u>	25x35	1.8 mil	Orange/Black	13-16 gal	200
<u>886</u>	38x47	1.8 mil	Orange/Black	44 gal	100
<u>888</u>	106x86	1.8 mil	Clear/Black	60+ gal	25
Print Label: No Print					
<u>5023</u>	38x47	1.8 mil	Clear/No Print	44 gal	100
<u>5034</u>	33x40	2 mil	Clear/Black	33 gal	200
<u>5207</u>	70x55	1.8 mil	Clear/No Print	60 gal	100

CAESA-LAB INC.

www.caesalabinc.com

info@caesalabinc.com

Tel (450) 441-5408