




Disposable plastic cassettes hold tissue specimens during the embedding process, as well as in a storage file. Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.



Item No.	Specification	Quantity (pcs)	Measurement (cm)
F1009W	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009Y	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009G	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009P	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009LB	Tissue Embedding Cassette	5000	48 X 32 X 44
F1009DB	Tissue Embedding Cassette	5000	48 X 32 X 44
F1010	Embedding Cassette	4000	48 X 32 X 32

Disposable tissue processing embedding cassettes resistant to the chemical action of histological solvents, anterior writing area at a 30° angle. These cassettes are used to process tissues or biopsies with a method called embedding process; this process is used to preserve tissues or biopsies with paraffin.



CAESA-LAB Inc.  
 4046 Rue de L'Orchidee  
 Saint-Bruno, QC J3V 6J9  
 Canada  
 Tel. (450)-441-5408  
 Fax (450)-441-9167  
[www.caesalabinc.com](http://www.caesalabinc.com)