

Champion S-33

Swing Out Rotor Bench-top Centrifuge

Centrifuge



Features

- Vortex ventilation system keeps sample temperature consistent , while maintaining sample integrity
- Wide mouth opening allows easy access to inner chamber
- Closed system inner chamber facilitates cleaning and maintenance
- Temperature controlled maintenance-free brushless motor
- Digital tachometer and RPM adjustment allows precise speed control
- Smooth acceleration and deceleration

Specifications

Item	Specification
Speed	300 - 3,300 rpm
Max Rcf	1732 xg
Timer	Digital, 0 – 30 min
Noise at Max Speed	≤65dB
Power	110V 60Hz or 220V 50Hz
Motor	Brushless
Weight	30 lb
Dimensions (LxWxH)	462 x 362 x 318 mm
Rotor	8 place duralumin Swing out rotor with 8x15ml PC buckets
Compatible Test Tubes	15ml test tubes w/ diameter less than 18mm