

SM10 Series

Zoom Stereo Microscope

Microscope



Features

- Zoom range of 0.75X-7.5X (Zoom ratio of 1:10)
- Continuous magnification from 7.5 – 75X on standard models
- Magnification dial located on side of microscope
- Fan-shaped base provides large working area
- Durable guided-rail stand secures microscope head in place and allows for one-hand operation
- Diopter adjustment on both eye tubes
- 4 configurations available

Specifications

| Model | SMZ10-B | SMZ10-BLL | SMZ10-T | SMZ10-TLL |
|-------------------------|------------------------|-------------------------------|-------------------------|-------------------------------|
| Specification | | | | |
| Head | 45° inclined Binocular | 45° inclined Binocular | 45° inclined Trinocular | 45° inclined Trinocular |
| Eyepiece | WF10X/20mm | | | |
| Interpupillary Distance | 50mm-75mm | | | |
| Zoom Range | 0.75X-7.5X | | | |
| Range of Magnification | 7.5-75X | | | |
| Working Distance | 102mm | | | |
| Visual Field | 26.7-2.7mm | | | |
| Illumination | No Illumination | Top & Bottom LED Illumination | No Illumination | Top & Bottom LED Illumination |
| Power | N/A | 90V-250V 50Hz/60Hz | N/A | 90V-250V 50Hz/60Hz |

Accessorize your SMZ10 series microscope



Fluorescent Ring Lights



LED Ring Light



Fiber Optics Illuminator



Boom Stand



CMOS Camera