

Tissue Cassettes in Dispenser Boxes



We offer two great cassettes – one is an alternate to the Sakura Uni-Cassette® which is ideal for automated cassette printers. The other is a high quality economy style for manual procedures.

Tissue Processing / Embedding Cassettes For Automated Processing

For use on most automated cassette printers. This cassette is a great alternative to the Sakura Uni-Cassette®. The one-piece design allows for the lids to be closed and reopened repeatedly without the lid getting detached from the base.

Packaging: 500/dispenser box, 2 boxes/case

35° angled writing area.

TISSUE CASSETTE		BIOPSY CASSETTE	
Item#	Color	Item#	Color
1094A	Aqua	1095B	Blue
1094B	Blue	1095G	Green
1094G	Green	1095N	Orange
1094H	Peach	1095P	Pink
1094L	Lilac	1095W	White
1094N	Orange	1095Y	Yellow
1094P	Pink		
1094T	Tan		
1094W	White		
1094Y	Yellow		


Tissue Processing / Embedding Cassettes For Manual Processing

Ideal for holding tissue specimens during the embedding process. The integral lid easily snaps apart and locks into the base of the cassette, which can be opened and closed repeatedly.

Packaging: 500/dispenser box, 2 boxes/case

30° angled writing area.

TISSUE CASSETTE		BIOPSY CASSETTE	
Item#	Color	Item#	Color
1090A	Aqua	1091B	Blue
1090B	Blue	1091G	Green
1090G	Green	1091W	White
1090P	Pink		
1090W	White		
1090Y	Yellow		



Formalin Filled Containers

We can offer you contract pricing so you can target large accounts and private label contracts!

Item#	Description	Unit
6518FL	20 mL (0.66 oz), 35 mm ID, 44 mm H	100/case
6520FL	40 mL (1.33 oz), 48 mm ID, 41 mm H	100/case
6522FL	60 mL (2 oz), 48 mm ID, 60 mm H	100/case
6524FL	90 mL (3 oz), wide-mouth, 48 mm ID, 76 mm H	100/case
6525FL	90 mL (3 oz), wide-mouth, 53 mm ID, 63 mm H	100/case
6527FL	120 mL (4 oz), wide-mouth, 53 mm ID, 78 mm H	100/case



Items #6518FL, 6520FL, 6522FL, 6524FL, 6525FL, 6527FL