

Inoculation Loops

INOCULATION LOOP - RIGID FEEL

These disposable inoculation loops and needles offer the most convenient and safest way of streaking culture media. Molded from a special medical grade of polystyrene, they offer a strong, rigid feel to provide comfortable, stable and uniform streaking. The loops are enclosed in peel-apart packages or individually wrapped. All loops are sterilized by gamma radiation. The 1 μ L and 10 μ L inoculation loops are available with the loop on one end and an integrated needle on the other end. A separate needle is therefore not required. We also offer a 10/1uL combination loop consisting of 1uL on one end and 10uL on the other end.

- Certificate of calibration
- Gamma Sterilized

Item #	Description	Unit
2851	1uL with needle, blue, STERILE, individually wrapped	500
2861	1uL with needle, blue, STERILE, 10/peel pack (50 packs per case)	500
2811	1uL with needle, blue, STERILE, 20/peel pack (25 packs per case)	500
2855	10uL with needle, yellow, STERILE, individually wrapped	500
2865	10uL with needle, yellow, STERILE, 10/peel pack (50 packs per case)	500
2815	10uL with needle, yellow, STERILE, 20/peel pack (25 packs per case)	500
2850	10/1uL, white, STERILE, individually wrapped	500
2860	10/1uL, white, STERILE, 10/peel pack (50 packs per case)	500
2810	10/1uL, white, STERILE, 20/peel pack (25 packs per case)	500



INOCULATION LOOP - FLEXIBLE FEEL

These inoculation loops are produced from medical grade polystyrene. Both the loop and needle ends are designed to be more bendable than our "rigid" loops. They are smooth and flexible to facilitate uniform streaking of culture media. These loops are ETO sterilized and packed 20 loops per peel-pack.

Item #	Description	Unit
2801	1uL with needle, natural, STERILE, 20/peel pack (50 packs per case)	1000
2805	10uL with needle, blue, STERILE, 20/peel pack (50 packs per case)	1000



QUADLOOP® INOCULATION LOOP - WITH SPREADING SPHERE

The new disposable Quadloop enables sampling and dilution streaking without the need to change loops. This is possible due to the unique sphere/handle arrangement which provides four sterile sampling and streaking surfaces simply by turning the handle! On the other end of the Quadloop® is a regular 1uL or 10uL loop for sampling fixed volumes, or a needle for stabbing.


- Gamma Sterilized

Item #	Description	Unit
2841	1uL with spreading sphere, green, STERILE, 20/peel pack (50 packs per case)	1000
2845	10uL with spreading sphere, blue, STERILE, 20/peel pack (50 packs per case)	1000
2813	Needle with spreading sphere, orange, STERILE, 20/peel pack (50 packs/case)	1000
2863	Needle with spreading sphere, orange, STERILE, 10/peel pack, 50 packs/case)	500



SPREADERS AND NEEDLES

Item #	Description	Unit
2910	L-shaped spreader, polystyrene, STERILE, individually wrapped	100
2911	T-shaped spreader, polystyrene, STERILE, individually wrapped	100
2813	Needle with spreading sphere, orange, STERILE, 20/peel pack (50 packs/case)	1000
2863	Needle with spreading sphere, orange, STERILE, 10/peel pack, 50 packs/case)	500



CAESA-LAB INC.
Tel: (450) 441-5408
www.caesalabinc.com