



# Microscope Slides

## MICROSCOPE SLIDES


Caesa-Lab premium quality slides are quickly replacing overpriced, name-brand slides in markets throughout the world. From the smooth finish of our chemically etched frosted slides, to our unique safety corners, Caesa-Lab slides will exceed your expectations. Produced from premium glass, Caesa-Lab slides are packaged in boxes of 72 and feature a moisture absorbent sheet in each box to prevent sticking. All slides are pre-cleaned and ready to use.



Standard Slide



Corners are slightly clipped to prevent finger and hand cuts



No sharp edges  
Perfect for blood smears

### Plain

Caesa-Lab's plain microscope slides are available in 3 styles: ground edge, ground edge with safety corners, and beveled edge with clipped corners.

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Packaging: 72/box, 20 boxes/case (10 gross)

Item #	Description	Unit
1301	Plain, 90° ground edges	1440
1321	Plain, 90° ground edges with safety corners	1440
1331	Plain, 45° beveled edges with clipped corners	1440




### Frosted

Caesa-Lab's frosted microscope slides are produced using a unique chemical etching process. This provides our slides with a smooth frosted finish that won't wear out slide marking pens. Slides are printed under the frosting area for orientation purposes. Finding the perfect slide to suit your needs is easy: choose from single sided or double sided frosted, and ground or beveled edges. We also offer our ground-edge slides with the following options: 90° corners (straight-edge) or safety corners. We are happy to help you select the perfect slide. Contact us at (201) 599-1400 for assistance.

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Packaging: 72/box, 20 boxes/case (10 gross)

Item #	Description	Unit
1304	Frosted 1 end, 1 side, 90° ground edges, 90° corners	1440
1324	Frosted 1 end, 1 side, 90° ground edges with safety corners	1440
1328	Frosted, 1 end, 2 sides, 90° ground edges with safety corners	1440
1334	Frosted 1 end, 1 side, 45° beveled edges with clipped corners	1440
1338	Frosted 1 end, 2 sides, 45° beveled edges with clipped corners	1440




# Microscope Slides

## **Frosted - Color**

Caesa's frosted microscope slides are offered in bright, attractive colors. They are available in three edge styles: ground, ground with safety corners and beveled. The ground edge slide offers 90° corners and is perfect for use in automatic cover slip machines. The ground edge slide with safety corners provide slightly clipped corners for comfortable handling and to prevent the risk of finger and hand cuts. The beveled edge slide is ideal for performing blood smears. The highly polished edges are specially angled to 45° to provide a secure grasp during the smearing process. All slides are printed under the frosting area for orientation purposes. They are pre-cleaned and ready to use.

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Frosted one end, one side

Packaging: 72/box, 20 boxes/case (10 gross)




## **Ground Edge**

Item #	Description	Unit
1304W	Frosted, 90° ground edge, WHITE	1440
1304B	Frosted, 90° ground edge, BLUE	1440
1304G	Frosted, 90° ground edge, GREEN	1440
1304N	Frosted, 90° ground edge, ORANGE	1440
1304P	Frosted, 90° ground edge, PINK	1440
1304Y	Frosted, 90° ground edge, YELLOW	1440




## **Ground Edge with Safety Corners**

1324W	Frosted, 90° ground edge with safety corners, WHITE	1440
1324A	Frosted, 90° ground edge with safety corners, AQUA	1440
1324B	Frosted, 90° ground edge with safety corners, BLUE	1440
1324D	Frosted, 90° ground edge with safety corners, GOLD	1440
1324G	Frosted, 90° ground edge with safety corners, GREEN	1440
1324L	Frosted, 90° ground edge with safety corners, LILAC	1440
1324N	Frosted, 90° ground edge with safety corners, ORANGE	1440
1324P	Frosted, 90° ground edge with safety corners, PINK	1440
1324T	Frosted, 90° ground edge with safety corners, TAN	1440
1324Y	Frosted, 90° ground edge with safety corners, YELLOW	1440



## **Beveled Edge with Clipped Corners**

1334W	Frosted, 45° beveled edge with clipped corners, WHITE	1440
1334B	Frosted, 45° beveled edge with clipped corners, BLUE	1440
1334G	Frosted, 45° beveled edge with clipped corners, GREEN	1440
1334N	Frosted, 45° beveled edge with clipped corners, ORANGE	1440
1334P	Frosted, 45° beveled edge with clipped corners, PINK	1440
1334Y	Frosted, 45° beveled edge with clipped corners, YELLOW	1440




# Microscope Slides

## **Charged Microscope Slides**

Caesa-Lab's charged microscope slides are offered in bright, attractive colors and provide the ideal surface for cell and tissue adhesion. These popular slides are mostly used in histology, cytology and pathology departments, where the convenience of an adhesive slide saves the user time and money. They are pre-cleaned and ready to use.

- Dimensions: 25 x 75 x 1mm
- Thickness: 1 - 1.2mm
- Frosted one end, one side

Packaging: 72/box, 20 boxes/case (10 gross)



## **Ground Edge**

Item #	Description	Unit
1358A	Charged microscope slide, Aqua	1440
1358B	Charged microscope slide, Blue	1440
1358D	Charged microscope slide, Gold	1440
1358G	Charged microscope slide, Green	1440
1358L	Charged microscope slide, Lilac	1440
1358N	Charged microscope slide, Orange	1440
1358P	Charged microscope slide, Pink	1440
1358T	Charged microscope slide, Tan	1440
1358W	Charged microscope slide, White	1440
1358Y	Charged microscope slide, Yellow	1440

## **Beveled Edge with Clipped Corners**

1354W	Charged microscope slide, beveled, White	1440
1354B	Charged microscope slide, beveled, Blue	1440
1354G	Charged microscope slide, beveled, Green	1440
1354N	Charged microscope slide, beveled, Orange	1440
1354P	Charged microscope slide, beveled, Pink	1440
1354Y	Charged microscope slide, beveled, Yellow	1440

## **Microscope Slides with Cavity**

Our 25 x 75 x 1mm cavity microscope slides are pre-cleaned and ready to use. They are offered with 90° ground edges and feature our popular safety corners. Two styles are available: single cavity or double cavity.

Item #	Description	Unit
1341	Microscope slide with single cavity, 72/box, 20 boxes/case	1440
1341-72	Microscope slide with single cavity, 72/box	72
1344	Microscope slide with double cavity, 72/box, 20 boxes/case	1440
1344-72	Microscope slide with double cavity, 72/box	72



# Microscope Slides

## COVER GLASS


Our cover glass is packaged: 1 oz per box, 10 boxes per carton (10 oz). Each flip top box contains a desiccant and is wrapped in aluminum to prevent sticking.

Available in three thicknesses:

No. 1: 0.14mm +/- 0.015mm

No. 1.5: 0.175mm +/- 0.015mm

No. 2: 0.21mm +/- 0.2mm



Item #	Description	Unit	Approximate Pcs per 1 oz box
1401-10	18 x 18mm, No. 1	10 oz	248
1401-20	18 x 18mm, No. 2	10 oz	168
1402-10	20 x 20mm, No. 1	10 oz	192
1402-20	20 x 20mm, No. 2	10 oz	154
1404-10	22 x 22mm, No. 1	10 oz	168
1404-15	22 x 22mm, No. 1.5	10 oz	148
1404-20	22 x 22mm, No. 2	10 oz	136
1407-10	22 x 32mm, No. 1	10 oz	114
1407-20	22 x 32mm, No. 2	10 oz	84
1411-10	22 x 40mm, No. 1	10 oz	86
1411-20	22 x 40mm, No. 2	10 oz	75
1414-10	22 x 50mm, No. 1	10 oz	70
1414-20	22 x 50mm, No. 2	10 oz	63
1418-10	22 x 60mm, No. 1	10 oz	60
1418-20	22 x 60mm, No. 2	10 oz	50
1405-10	24 x 24mm, No. 1	10 oz	141
1405-20	24 x 24mm, No. 2	10 oz	102
1408-10	24 x 32mm, No. 1	10 oz	104
1408-20	24 x 32mm, No. 2	10 oz	78
1412-10	24 x 40mm, No. 1	10 oz	82
1412-20	24 x 40mm, No. 2	10 oz	71
1415-10	24 x 50mm, No. 1	10 oz	65
1415-15	24 x 50mm, No. 1.5	10 oz	56
1415-20	24 x 50mm, No. 2	10 oz	56
1419-10	24 x 60mm, No. 1	10 oz	55
1419-20	24 x 60mm, No. 2	10 oz	46

# Microscope Slides

## **SLIDE MAILERS**

Caesa's sturdy cardboard and plastic slide mailers are designed for the safe, easy handling of microscope slides. The 1, 2 and 3 place slide mailers securely hold 25 x 75mm and 1" x 3" slides with or without a cover-glass. They feature an inner finger groove for effortless slide removal along with a cover that opens for viewing. The flat surface area is ideal for labeling or writing any necessary information. They are light-weight and reusable to minimize cost.

The 20-place cardboard slide holder contains a numerical index that corresponds to the inner slide compartments. The inside flaps are recessed so that there is no contact or contamination between the slides and the cover when closed. With a 360° hinge, the cover folds flat underneath to save space when the holder is in use.





## **Cardboard Mailers**

Item #	Description	Unit
513001	Cardboard mailer for 1 slide, 100/box, 1000/case	1000
513011	Cardboard mailer for 1 slide, 100/case	100
513002	Cardboard mailer for 2 slides, 50/box, 1000/case	1000
513022	Cardboard mailer for 2 slides, 50/box, 100/case	100
513003	Cardboard mailer for 3 slides, 50/box, 500/case	500
513013	Cardboard mailer for 3 slides, 50/case	50
513004	Cardboard mailer for 4 slides, 50/box, 200/case	200
513014	Cardboard mailer for 4 slides, 50/case	50
513025	Cardboard mailer for 20 slides, 10/case	10
513026	Cardboard mailer for 20 slides, Each	1




## **Plastic Mailers - Extra Sturdy**

Item #	Description	Unit
513031	Plastic mailer for 1 slide, 100/box, 1000/case	1000
513032	Plastic mailer for 1 slide, 100/box	100
513041	Plastic mailer for 2 slides, 100/box, 1000/case	1000
513042	Plastic mailer for 2 slides, 100/box	100
513045	Plastic mailer for 3 slides, 100/box, 1000/case	1000
513046	Plastic mailer for 3 slides, 100/box	100



## **Plastic Mailers with Flip Top Lid**

Item #	Description	Unit
513051	Plastic mailer for 2 slides, natural, 100/box	100
513052	Plastic mailer for 2 slides, natural, 10 x 100/case	1000
513051B	Plastic mailer for 2 slides, blue, 100/box	100
513052B	Plastic mailer for 2 slides, blue, 10 x 100/case	1000
513061	Plastic mailer for 5 slides, natural, 100/box	100
513062	Plastic mailer for 5 slides, natural, 10 x 100/case	1000



# Microscope Slides

## **MICROSCOPE SLIDE STORAGE BOXES**

Caesa's storage boxes are designed for the safe long term storage of microscope slides. They are produced from durable ABS plastic and offer a hinged lid for easy opening and closing. An inventory sheet is located on the inside cover for convenient slide identification and organization. The base of the box is cork lined for added protection. Available in four sizes, these storage boxes accommodate 25 x 75mm and 1 x 3" microscope slides.



<b>Item #</b>	<b>Description</b>	
513072A	Storage box for 12 microscope slides	1
513075A	Storage box for 25 microscope slides	1
513077A	Storage box for 50 microscope slides	1
513079A	Storage box for 100 microscope slides	1

## **SLIDE DISPENSER**

The Caesa slide dispenser is designed for use with 25 x 75mm and 3 x 1" microscope slides. To use, simply insert one full box of 72 slides into the holder. Turn the dial to dispense one slide at a time. This prevents slides from getting dirty, dusty and scratched.

<b>Item #</b>	<b>Description</b>	<b>Unit</b>
513000	Microscope slide dispenser for 72 slides	1

