


# Microtest Plates

## 96-WELL MICROTEST PLATES


Three types of 96-well microtest plates are offered for general (non-tissue culture) applications including microbiology and serology. Available in "U" (round bottom), "V" and flat bottom shapes, they are available individually wrapped to protect them from air contaminants. Made from high optic transmission grade polystyrene, the microtest plates are arranged in a standard 12 x 8 configuration and have alpha-numeric position markings. Optional lids are also available. Both plates and lids may be purchased sterile or non-sterile.



- Individually wrapped to protect from air contaminants.
- Available in "U", "V" and Flat bottoms


U-Bottom	V-Bottom	Flat Bottom
Each well is 300uL	Each well is 250uL	Each Well is 375uL

Item #	Description	Unit
120030	U-bottom, individually wrapped, non-sterile	50
120038	U-bottom, individually wrapped, STERILE	50
120130	V-bottom, individually wrapped, non-sterile	50
120138	V-bottom, individually wrapped, STERILE	50
120330	Flat-bottom, individually wrapped, non-sterile	50
120338	Flat-bottom, individually wrapped, STERILE	50
129930	Lid, individually wrapped, non-sterile	150
129938	Lid, individually wrapped, STERILE	150



## SEALING FILM

Designed for the sealing of microplates, our adhesive sealing films eliminate well-to-well contamination and crossover in sensitive ELISA/EIA assays and PCR applications where minimization of evaporation and vapor loss is critical. Additionally, our films can be used to seal tissue culture plates for short term storage and incubation, transport, and the containment of biohazardous solutions. An economical alternative to plate covers, adhesive sealing films provide a complete seal, whereas plate lids act only as covers.



Item #	Description	Unit
129940	EZ-Seal sealing film, non-sterile	100
129948	EZ-Seal sealing film, STERILE	100
129950	Breathable sealing film, non-sterile	100
129958	Breathable sealing film STERILE	100
129980	Tuff-Seal extra strong film, non-sterile	100
129990	Aluminum film, non-sterile	100
129946	Mini-Strips sealing film, non-sterile	200