


# Racks

## RACKS FOR TEST TUBES

Our well-designed test tube racks are made of high density polypropylene, and are autoclavable and non-floating in a water bath. Sizes range from 13mm to 30mm and, therefore, accommodate all routine tubes. The 30mm size rack is especially suitable for 50mL centrifuge tubes. Alphanumeric markings are molded into the racks to provide easy test tube identification. They are available in four colors for color coding of test samples, or simply for aesthetics in the modern laboratory. The racks are shipped flat in order to reduce shelf space and then simply snapped together for use. They are stackable for convenient storage.

- Produced from polypropylene
- Autoclavable at 12 1°C for 20 minutes
- Non-floating in a water bath.
- Can be used in freezers without deforming or becoming brittle



### **For 13mm Tubes**

90 holes, 6 x 15 format, rack dimensions: 246 x 105 x 64mm

Item #	Description	Unit
456403	13mm size, 90-place, white	Each
456404	13mm size, 90-place, blue	Each
456406	13mm size, 90-place, yellow	Each
456410	13mm size, 90-place, red	Each

### **For 16mm Tubes**

60 holes, 5 x 12 format, rack dimensions: 246 x 105 x 72mm

#### **For 25mm Tubes**

Item #	Description	Unit
456503	16mm size, 60-place, white	Each
456504	16mm size, 60-place, blue	Each
456506	16mm size, 60-place, yellow	Each
456510	16mm size, 60-place, red	Each

### **For 20mm Tubes**

40 holes, 4 x 10 format, rack dimensions: 246 x 105 x 72mm

Item #	Description	Unit
456603	20mm size, 40-place, white	Each
456604	20mm size, 40-place, blue	Each
456606	20mm size, 40-place, yellow	Each
456610	20mm size, 40-place, red	Each

## RACKS FOR MICROCENTRIFUGE TUBES

Caesa's sturdy microcentrifuge tube racks are made of high density polypropylene and hold up to 100 microcentrifuge tubes. They are autoclavable and are non-floating in a water bath. Alphanumeric grid marks are molded into the racks to identify tube positions. The racks are supplied flat to reduce shelf space and transportation costs. Racks are stackable and have spaces for labels on each side (for bar coding, etc...)

- 100 holes, rack dimensions: 263 x 109 x 45mm

Item #	Description	Unit
456303	Microcentrifuge Tube Rack, 100-place, white	Each

40 holes, 4 x 10 format, rack dimensions: 295 x 125 x 85mm

Item #	Description	Unit
456703	25mm size, 40-place, white	Each
456704	25mm size, 40-place, blue	Each
456706	25mm size, 40-place, yellow	Each
456710	25mm size, 40-place, red	Each

### **For 30mm Tubes**

24 holes, 3 x 8 format, rack dimensions: 300 x 112 x 85mm

Item #	Description	Unit
456803	30mm size, 24-place, white	Each
456804	30mm size, 24-place, blue	Each
456806	30mm size, 24-place, yellow	Each
456810	30mm size, 24-place, red	Each

