

SPECIMEN TRANSPORT BAGS

SPECIMEN TRANSPORT BAGS - ZIPPER CLOSURE

Our specimen transport bags are designed for the safe transportation of specimens and documentation. Each bag is printed with the biohazard symbol and offers checkmark boxes to indicate special handling instructions.

- 2 mil gauge polyethylene
- Attached document compartment
- Zipper closure allows repeated access to specimen
- Imprinted with biohazard symbol

Item #	Description	Unit
4918	6" x 9" specimen transport bag, zipper closure	1000
4920	6" x 9" specimen transport bag with tear-zone, zipper closure	1000
4925	8" x 10" specimen transport bag with tear-zone, zipper closure	1000
4928	12" x 15" specimen transport bag with tear-zone, zipper closure	1000
4950A	6 x 9" STAT specimen transport bag, red tint, zipper closure	1000
4950	6 x 9" STAT specimen transport bag with tear-zone, red tint, zipper closure	1000



#4918



#4920

Custom printed bags are available. Please contact our sales department for information.

SPECIMEN TRANSPORT BAGS - GLUE CLOSURE

The special glue closure of our specimen transport bags offer a liquid-tite tamper evident seal. Simply pull off the liner to expose the glue adhesive and press to close. Item #4923 features a 3 x5" liquid absorbent pad that quickly absorbs aqueous solutions or bodily fluids providing for safe transportation of specimens.

- 1.5 mil gauge polyethylene
- Attached document compartment
- Glue closure tamper evident seal
- Imprinted with biohazard symbol
- Perfect for the transport of drug testing products

Item #	Description	Unit
4921	6" x 9" specimen transport bag, glue closure	1000
4923	6" x 9" specimen transport bag, glue closure, with absorbent pad.	500



BITRAN SERIES S STORAGE BAGS

Bitran Series S bags are designed for long-term handling and storage of specimens. The double track zipper is easy to close and provides an airtight seal. Produced from 3 mil Saranex, these bags provide exceptional vapor and odor control.

- For long-term storage of specimens
- Appropriate for use as a primary container

Item #	Description	Unit
4740	2" x 4", 4 x 250/case	1000
4741	3" x 6", 4 x 250/case	1000
4742	6" x 6", 4 x 250/case	1000
4743	7" x 8", 4 x 250/case	1000
4744	9" x 12", 4 x 250/case	1000
4745	12" x 12", 4 x 250/case	1000
4746	16" x 16", 100/case	100

