

# 30ml Collectors

for specimen collection and transportation



multi-purpose collector

slide holder

fecal collection

pinworm collection

s a f e t y , q u a l i t y

& r e l i a b i l i t y

# 30ml Collectors

## the vial

Each vial is manufactured using medical grade plastic which ensures the consistency and reliability of the product. Each style of container is leakproof and travels well during harsh transport or within a pneumatic tube system. The Fecal, Pinworm and Universal collectors have graduations on the vial wall for easy measurement. The full circle thread design of this screw cap container ensures that no cross threading will occur.

- 30ml clarified polypropylene
  - break-resistant
  - graduated
  - skirted conical bottom
  - leakproof
  - labeled
- easy visual analysis of sample  
safe during transport and sampling  
easy quantification of sample  
convenient, free standing  
protection against contamination  
ample space for patient information

## the cap

- “o-ring” cap design
  - ridges on cap
- leakproof closure  
allows secure grip

## the fecal collector

Features a combination spoon and fork device to allow easy collection of any form of fecal matter. The vial has easy to read graduations which assist in the collection of the proper sample size.

## the pinworm collector

Features a clear polystyrene paddle, lightly adhesive on one side allowing easy collection of ova.

## the universal collector

This versatile, freestanding, skirted conical bottom vial can be used for multiple purposes including fecal, urine and sputum collection.

## the slide holder

Ideal for transporting up to four standard size slides with or without fixative. Slides do not touch and remain protected from breakage within each groove.  
This unlabeled product is suitable for use as a mailer.

code	description
V302	30ml clarified, graduated, polypropylene collector, white “o-ring” cap, labeled, sterile, 200/bag, 1000/case
V302-F	30ml clarified, graduated, polypropylene vial with fecal spork insert cap, labeled, 200/bag, 1000/case
V302-PIN	30ml clarified, graduated, polypropylene vial with pinworm paddle insert cap, labeled, bulk, 150/case
V302-SH	30ml clarified, polypropylene vial with groove inserts for four slides, red “o-ring” cap, no label, 250/bag, 500/case

## and more!

Ask us about other product options, such as our full line of 5ml and 10 ml polypropylene vials.

## CAESA-LAB

Telephone: (450) 441-5408

Fax (450) 441-9167

[www.caesalabinc.com](http://www.caesalabinc.com)

[info@caesalabinc.com](mailto:info@caesalabinc.com)