

Multi-Purpose Containers

NEW from Starplex

Innovation in specimen collection



Translucent containers for the storage of pathology, industrial, and environmental samples. These containers are made from high quality plastic. The matching caps snap on tight for a secure fit.

Available in sizes 4oz to 172oz (120ml to 5.1L). Bulk-packed.

- ◆ Economical
- ◆ Versatile
- ◆ Wide range of sizes
- ◆ Flexible & Durable

| Part. No. | Volume | Description | Units/Case |
|-----------|--------------|---------------------------------------------|------------|
| PA0610-T | 4oz / 120ml | Non-sterile container, snap lid | 250 |
| PA0620-T | 8oz / 240ml | Non-sterile container, snap lid | 250 |
| PA0640-T | 16oz / 480ml | Non-sterile container, snap lid | 100 |
| PA0670-T | 32oz / 950ml | Non-sterile container, snap lid | 100 |
| PA0690-T | 64oz / 1.9L | Non-sterile container, snap lid | 50 |
| PA0695-T | 85oz / 2.5L | Non-sterile container, snap lid | 50 |
| PA0700-T | 160oz / 4.8L | Non-sterile container with handle, snap lid | 25 |
| PA0720-T | 172oz / 5.1L | Non-sterile container with handle, snap lid | 25 |

CAESA-LAB INC.

TELEPHONE : (450) 441-5408

FAX : (450) 441-9167

E-MAIL : info@caesalabinc.com

www.caesalabinc.com