

TRANSFER PIPETS


P200

The DROPETTE®

DISPOSABLE TRANSFER PIPETS

These unbreakable all-in-one pipets eliminate the hazard of broken glass and exposure to infectious materials. Put an end to matching rubber bulbs with glass pipets. Molded from see-through low density polyethylene. Inert to biological fluids and most acids. The low-affinity surface reduces the loss of cells and valuable proteins due to binding. Can be sealed and refrigerated. They work well whenever there is a need for quick, safe transfer of fluids. Temperature resistant down to -196 °C. Can be gaz (EtO) sterilized.

Choose between 12 very popular models, available in several sizes, tip designs and lengths, in sterile or non sterile packaging. Seven models provide graduations.



- Will not shatter or break
- Can be used in liquid nitrogen
- Unbreakable
- Non toxic and inert
- No bulb to insert or remove
- Uniform drop size

Cat. #	Sterile	Graduated	Length	Capacity	Bulb draw	Drops per ml.	Inner pack
P200-10 P200-101S P200-1020S	• •		6.3 cm	1.2 ml	0.9 ml	25	Loose 1 20
P200-14 P200-141S P200-1420S	• •	•	11.4 cm	1.5 ml	0.8 ml	23	Loose 1 20
P200-20 P200-201S P200-2020S	• •		8.7 cm	1.7 ml	0.9 ml	25	Loose 1 20
P200-30 P200-301S P200-3020S	• •	•	13.8 cm	3 ml	2.1 ml	25	Loose 1 20
P200-44 P200-441S P200-445S P200-4410S P200-4420S	• • • •		15.5 cm	4 ml	3.1 ml	25	Loose 1 5 10 20
P200-52 P200-521S P200-525S P200-5210S P200-5220S	• • • •	•	15 cm	5 ml	3.1 ml	25	Loose 1 5 10 20
P200-56 P200-561S P200-565S P200-5610S P200-5620S	• • • •	•	15.5 cm	5 ml	1.9 ml	25	Loose 1 5 10 20
P200-58 P200-581S P200-5820S	• •	•	15.6 cm	5 ml	1.8 ml	25	Loose 1 20
P200-58V P200-58V1S P200-58V20S	• •	•	14.5 cm	5 ml	1.8 ml	25	Loose 1 20
P200-60 P200-601S P200-605S P200-6010S P200-6020S	• • • •		22.5 cm	6 ml	2.3 ml	25	Loose 1 5 10 20
P200-72 P200-721S P200-725S P200-7210S P200-7220S	• • • •	•	15.5 cm	7 ml	3.2 ml	22	Loose 1 5 10 20
P200-82 P200-821S P200-8220S	• •		15 cm	8 ml	4.6 ml	22	Loose 1 20

Packaging:

All non sterile pipets are in boxes of 500 and cases of 5000.

Exception: P200-82 are in boxes of 400 and cases of 4000.

Most sterile pipets are in boxes of 400 and cases of 4000. Please call Customer Service for more details.

