

30ml Collectors

for specimen collection and transportation



multi-purpose collector

slide holder

fecal collection

pinworm collection

s a f e t y , q u a l i t y
& r e l i a b i l i t y

Registered to ISO 9001:2000

ST RPLEX[®]
scientific inc.

30ml Collectors

the vial

Each vial is manufactured using medical grade plastic which ensures the consistency and reliability of the product. Each style of container is leakproof and travels well during harsh transport or within a pneumatic tube system. The Fecal, Pinworm and Universal collectors have graduations on the vial wall for easy measurement. The full circle thread design of this screw cap container ensures that no cross threading will occur.

- 30ml clarified polypropylene
 - break-resistant
 - graduated
 - skirted conical bottom
 - leakproof
 - labeled
- easy visual analysis of sample
safe during transport and sampling
easy quantification of sample
convenient, free standing
protection against contamination
ample space for patient information

the cap

- “o-ring” cap design
 - ridges on cap
- leakproof closure
allows secure grip

the fecal collector

Features a combination spoon and fork device to allow easy collection of any form of fecal matter. The vial has easy to read graduations which assist in the collection of the proper sample size.

the pinworm collector

Features a clear polystyrene paddle, lightly adhesive on one side allowing easy collection of ova.

the universal collector

This versatile, freestanding, skirted conical bottom vial can be used for multiple purposes including fecal, urine and sputum collection.

the slide holder

Ideal for transporting up to four standard size slides with or without fixative. Slides do not touch and remain protected from breakage within each groove.
This unlabeled product is suitable for use as a mailer.

starplex code	description
V302	30ml clarified, graduated, polypropylene collector, white “o-ring” cap, labeled, sterile, 200/bag, 1000/case
V302-F	30ml clarified, graduated, polypropylene vial with fecal spork insert cap, labeled, 200/bag, 1000/case
V302-PIN	30ml clarified, graduated, polypropylene vial with pinworm paddle insert cap, labeled, bulk, 150/case
V302-SH	30ml clarified, polypropylene vial with groove inserts for four slides, red “o-ring” cap, no label, 250/bag, 500/case

and more!

Ask us about other product options, such as our full line of 5ml and 10 ml polypropylene vials.

vision statement

Starplex Scientific Inc. is dedicated to the development of innovative, premium quality plastic disposables; providing reliable specimen transportation for many facets of the health care industry. Focused on and responsive to the needs of the customer, Starplex strives to anticipate the evolving market requirements and work co-operatively with customers to implement unique, effective solutions.

record of reliability

Our origins date back to the 70's, and since then we have become a prominent supplier to the medical industry. Our sales and distribution coverage includes every state and province of North America, supported by a state-of-the-art manufacturing plant in Canada. We are privileged to supply literally thousands of hospitals, laboratories and health care facilities with Starplex products.



50 A Steinway Blvd., Etobicoke, Ontario, Canada M9W 6Y3
tel: (416) 674-7474 / fax: (416) 674-6067 / toll free: 1-800-665-0954
email: info@starplexscientific.com / web: www.starplexscientific.com