


# Histoplex™

Tissue cassettes with micromesh chamber



Effective for the embedding, processing and storage of both standard tissue specimens and tissue fragments. The innovative micromesh chamber securely retains fragments leaving ample room available for additional larger specimens. Full range of colors available for identification.

s a f e t y , q u a l i t y

& r e l i a b i l i t y

**STARPLEX™**  
scientific inc.

# Histoplex™

## Features:

Durable acetyl polymer resists the chemical action of histological agents

Flow through slots maximize fluid exchange


Fine micromesh chamber, .23mm, ideal for tissue fragments, and effective for fluid exchange, no need for the additional use of filters

Textured, 45° angled writing surface, suitable for labeling devices

Clam-shell design can be snapped closed repeatedly and reopens with ease

Non-cytotoxic, range of 12 non-metallic colors

Conveniently packaged, 500 cassettes per dispenser box, 3 boxes per case



TC-WHITE

## Product Codes for Available Colors:

TC-AQUA	TC-BLUE	TC-YELLOW	TC-PINK	TC-ORANGE	TC-GREEN
					
TC-TAN	TC-RED	TC-LILAC	TC-GREY	TC-DRKGREEN	TC-PEACH
					

\*White also available

### vision statement

Starplex Scientific Inc. is dedicated to the development of innovative, premium quality plastic disposables; providing reliable specimen transportation for many facets of the health care industry. Focused on and responsive to the needs of the customer, Starplex strives to anticipate the evolving market requirements and work co-operatively with customers to implement unique, effective solutions.

### record of reliability

Our origins date back to the 70's, and since then we have become a prominent supplier to the medical industry. Our sales and distribution coverage includes every state and province of North America, supported by a state-of-the-art manufacturing plant in Canada. We are privileged to supply literally thousands of hospitals, laboratories and health care facilities with Starplex products.



50 Steinway Blvd., Etobicoke, Ontario, Canada M9W 6Y3  
 tel: (416) 674-7474 / fax: (416) 674-6067 / toll free: 1-800-665-0954  
 email: [info@starplexscientific.com](mailto:info@starplexscientific.com) / web: [www.starplexscientific.com](http://www.starplexscientific.com)